


# CONNECT

## PROVIDER PORTAL

## ERAD PACS

### EVO IMAGE VIEWER

GET REPORTS | VIEW IMAGES | DOWNLOAD IMAGES | ORDER EXAMS | MAKE APPOINTMENTS



**Main Menu**

**Window Level Toolbar**

**Toolbox Toolbar**

**Orientation Toolbar**

**Hanging Protocol Toolbar**

**Thumbnail Panel**

**Annotation Toolbar**

**Image Information Icon**

Key images appear in the thumbnail panel and have a key image icon on the top.







**Series Description**

**Visible in Frame**





**Unviewed Stack/Image**

**Series Number**










## MOUSE CONTROL

		<b>Scroll</b>	Press Left click and hold on lateral border of image, move mouse up and down to scroll through image series.
		<b>Window/Level</b>	Press Right click and hold on image, move mouse up and down, or left and right to increase or decrease window level.
		<b>Zoom</b>	Press Right and Left click together and hold on image, move mouse up to zoom in and down to zoom out. Press F1 to see full Page View.











## ORIENTATION TOOLBAR

	<b>Flip Horizontal</b>	Flip image(s) in selected frame horizontally.
	<b>Flip Vertical</b>	Flip image(s) in selected frame vertically.
	<b>Rotate Right</b>	Rotate image(s) in selected frame to the right 90°.
	<b>Rotate Left</b>	Rotate image(s) in selected frame to the left 90°.



## ANNOTATION TOOLBAR

	<b>Default Mode</b>	General cursor mode, used to select an image or object within the display.
	<b>Measure Distance</b>	Perform linear measurements.
	<b>Measure Angle</b>	Perform angular measurements.
	<b>Hounsfield Value</b>	Calculate and display the average pixel value within the defined.
	<b>Region of Interest</b>	Calculate the area within the defined region of interest.
	<b>Free Region</b>	Calculate the area within the defined region of interest.
	<b>Pointer</b>	Add a point graphic and optional text string.
	<b>Text</b>	Define a text string and place it on the image.
	<b>Spine Labeling</b>	Display labels identifying vertebral disks and interspaces.






## TOOLBOX TOOLBAR

	<b>Show Image Information</b>	Show/hide the demographic and image overlay information.
	<b>Patient Folder</b>	Access Patient Folder to View Reports and Prior Images for comparison.
	<b>Link Mode</b>	Link one or more image frames together.
	<b>Magic X</b>	Locate a defined point in one image on all other images that share the same frame of reference.
	<b>Magic Glass</b>	Display a floating magnifying glass for zooming in on the image data behind it.
	<b>Fit to Window</b>	Fit the selected image to the size of the frame.
	<b>Invert Grayscale</b>	Invert the slope of the applied window level.
	<b>Reset Orientation</b>	Reset the image(s) in the selected frame to the original orientation.
	<b>Reset Window/Level</b>	Reset the image(s) in the selected frame.
	<b>Localizer Lines</b>	The intersection of a referenced image on other views. Show All or Hide all localizer lines.

## HANGING PROTOCOL TOOLBAR

		Display applied hanging protocol. Expands hanging protocol selection list.
	<b>Create Layout</b>	Save the current layout as a hanging protocol template.
	<b>Layout Edit Mode</b>	Enable mode to customize layout frames.
	<b>Lock Grid Lines</b>	Enable or disable the ability to manually resize grid frames.
	<b>Tile Matrix</b>	Apply a tile pattern to a selected frame controlling the number of images displayed in the frame.
	<b>Grids Matrix</b>	Apply selected grid fill mode.
	<b>Grid fill Mode List</b>	Display selected grid fill mode and expands grid fill mode list.
	<b>Empty Cells</b>	Draw empty grid frames.
	<b>Reused Cells</b>	Draw grid frames and populate with currently displayed images.
	<b>Selected Study</b>	Draw grid frames and populate with selected study.
	<b>All Study</b>	Draw grid frames and populate with all studies loaded in viewer.



## WINDOW/LEVEL TOOLBAR

	Window width of the selected image frame.
	Window level of the selected image frame.
	Menu of preset window/level settings.
	Save the current W/L settings as a modality specific preset.
	The image enhancement list provides access to preset image sharpening, smoothing and edge detection.

## CT PRESET WINDOW/LEVEL HOTKEYS

PRESET	TOP ROW KEYBOARD
Initial	0
Soft	1
Liver	2
Lung	3
Brain	4
Subdural	5
Bone	6
Actual	8
Auto	9

## SESSION

	<b>Save</b>	Save Image or Series.
	<b>Print</b>	Print Image or Series.